

# GEN-SAS<sup>TM</sup>

## *Sperm Analyzer*

Actual screen magnification  
**400X, 800X, 1600X**

### **Easy**

Easy to operate

Rapid Semen Analysis

### **Excellent**

Three in one

Useful for multi-purposes

### **Economic**

SA without Makler Chamber



**IVF Gen**

*Contribution 2 next Generation*

IVF GEN Australia- Phone : +61 413 665 905

Email : [info@ivfgen.com](mailto:info@ivfgen.com), [ivfgen@gmail.com](mailto:ivfgen@gmail.com), Website : [www.ivfgen.com](http://www.ivfgen.com)

# GEN-SAS™

## Sperm Analyzer



Microscope



Makler Chamber



## Specifications

Product	GEN-SAS™
Dimensions(WxDxH)	150x200x330mm(Monitor not included)
Monitor	12 inch TFT LCD Monitor(With grid)
Camera	1/3" CCD Camera
Magnification	40X
Objective lens	X10, X20, X40
Plan stage	115(W) x 115(D)mm
Mechanical	Up and Down, X/Y axis, L type mechanical
Illumination	White LED
Power	DC 12V
Accessories	AV Cable, Power adapter

\* To improve the performance of equipment specifications may change without notice.



**IVF Gen**

*Contribution 2 next Generation*

**IVF GEN Australia- Phone : +61 413 665 905**  
**Email : info@ivfgen.com, ivfgen@gmail.com**  
**Website : www.ivfgen.com**