


IVF Gen- Micropipettes Product Catalogue



IVF Gen

Contribution 2 next Generation

IVF Gen ICSI Pipettes (GEN MIC) 10 Pcs/Box

Product Code	ID(μ m)	Tip Angle	Bevel Angle	Toe Length in mm (4 lengths available)
GEN MIC4525	4.5	25°	22° - 27°	550/650/750/1000
GEN MIC4525	4.5	30°	22° - 27°	550/650/750/1000
GEN MIC4525	4.5	35°	22° - 27°	550/650/750/1000
GEN MIC4525	5.0	25°	28° - 30°	550/650/750/1000
GEN MIC4525	5.0	30°	28° - 30°	550/650/750/1000
GEN MIC4525	5.0	35°	28° - 30°	550/650/750/1000

IVF Gen Holding Pipettes (GEN MHP) 10 Pcs/Box

Product Code	OD(μ m)	ID(μ m)	Tip Angle	Toe Length in mm (4 lengths available)
GEN MHP9030	90	15-20	30°	550/650/750/1000
GEN MHP9035	90	15-20	35°	550/650/750/1000
GEN MHP10030	100	15-20	30°	550/650/750/1000
GEN MHP10035	100	15-20	35°	550/650/750/1000
GEN MHP12030	120	15-20	30°	550/650/750/1000
GEN MHP12035	120	15-20	35°	550/650/750/1000
GEN MHP120B30	120	25-30	30°	550/650/750/1000
GEN MHP13035	130	15-20	35°	550/650/750/1000

IVF Gen Blastomere Biopsy Pipettes (GEN MBB) 10 Pcs/Box

Product Code	ID(μm)	Tip Angle	Toe Length in mm (4 lengths available)
GEN MBB182230	18-22	30°	550/650/750/1000
GEN MBB232730	23-27	30°	550/650/750/1000
GEN MBB283230	28-32	30°	550/650/750/1000
GEN MBB333730	33-37	30°	550/650/750/1000
GEN MBB384230	38-42	30°	550/650/750/1000
GEN MBB1822X30	18-22	30°	550/650/750/1000
GEN MBB2327X30	23-27	30°	550/650/750/1000
GEN MBB2832X30	28-32	30°	550/650/750/1000
GEN MBB3337X30	33-37	30°	550/650/750/1000
GEN MBB3842X30	38-42	30°	550/650/750/1000

IVF Gen Zona Drilling Pipettes (GEN MZD) 10 Pcs/Box

Product Code	ID(μ m)	Tip Angle
GEN MZD1015	8-10	15
GEN MZD1020	8-10	20
GEN MZD1025	8-10	25
GEN MZD1030	8-10	30
GEN MZD1035	8-10	35
GEN MZD1040	8-10	40
GEN MZD1045	8-10	45
GEN MZD1200	10-12	0
GEN MZD1205	10-12	5
GEN MZD1210	10-12	10
GEN MZD1215	10-12	15
GEN MZD1220	10-12	20
GEN MZD1225	10-12	25
GEN MZD1230	10-12	30
GEN MZD1235	10-12	35

IVF Gen Denuding Pipettes (GEN MDP) 10 Pcs/Box

PRODUCT CODE	ID(μm)
GEN MDP70	70 \pm 5 μm
GEN MDP100	100 \pm 5 μm
GEN MDP120	120 \pm 5 μm
GEN MDP130	130 \pm 5 μm
GEN MDP140	140 \pm 5 μm
GEN MDP150	150 \pm 5 μm
GEN MDP160	160 \pm 5 μm
GEN MDP170	170 \pm 5 μm
GEN MDP180P	180 \pm 5 μm
GEN MDP190P	190 \pm 5 μm
GEN MDP200P	200 \pm 5 μm
GEN MDP210P	210 \pm 5 μm
GEN MDP220P	220 \pm 5 μm
GEN MDP230P	230 \pm 5 μm
GEN MDP240P	240 \pm 5 μm
GEN MDP250P	250 \pm 5 μm
GEN MDP290P	290 \pm 5 μm
GEN MDP320P	320 \pm 5 μm
GEN MDP350P	350 \pm 5 μm



IVF Gen

Contribution 2 next Generation

IVF GEN Australia - Phone : +61 413 665 905

Email : info@ivfgen.com, ivfgen@gmail.com, Website : www.ivfgen.com